REDCap – Enable an Instrument as a Survey

• Step 1: On the Project Setup page and under "Main project settings," select the Enable button next to "Use surveys in this project?"

A P	Project Home	j≡ Project Setup	Other Functionality	Project Revision History	
Project status: 🖋 Development Completed steps 0 of 7					
	Mai	n project settings			
Enable 😔 Use surveys in this			in this project? ?	VIDEO: How to create and managed	ge a survey
Not st	Enat	Enable 😔 Use longitudinal data collection with defined events? ?			
l'm d	lone! Mo	dify project title, purp	oose, etc.		

• Step 2: Create your data collection instrument. Once you have created your instrument, go to the Online Designer and select the Enable button next to your instrument name, under the "Enabled a survey" column header.

The Online Designer will allow you to make project modifications to fields and data collection instruments very easily using only your web browser. NOTE: While in development status, all field changes will take effect immediately in real time.

Data Collection Survey options: Instruments Image: Survey Queue Image: Survey Survey Notifications		Survey options: I≡ Survey Queue Survey Login Survey Notifications	Ad Q	d new Create Import	instrumen a new ins a new ins a new ins	t: trument from scratch strument from the off ent ZIP file from anoth	n ficial <u>REDCap Shared Library</u> (e) her project/user or <u>external libraries</u> (e)
	Instrument name		Fields	View PDF	Enabled a survey	Instrument actions	Survey-related options
My First Instrument		4	7	Enable	Choose action 🗢		
	survey 2		1	7	Enable	Choose action 🗢	

• This will bring you into the Survey Settings, where you can give your survey a title, add survey instructions, edit the survey design options, and add other survey customizations.

A D		
T Pro	ect Hom	e

🚍 Project Setup 🛛 🍮 Online Designer

💿 Set Up My Survey

You may utilize any data collection instrument as a survey by enabling the instrument on this page. Data may still be collected on the instrument via normal data entry on the form, but enabling it as a survey allows for the additional option of collecting data via survey by respondents.

Set up my survey for data collection instrument "Demographics" Cancel				
Basic Survey Options:				
📎 Survey Title	My First Instrument			
Survey Instructions (Displayed at top of survey after title)	Please complete the survey below. Thank you!			
	B I U ANG ▲ · 型 · E E E E E - E E E E E C @ @ Mmm X & C @ @ M = E E E F C W @ Ø Mmm Y How to use Piping here			
Survey Design Options:				
Logo (Optional: display an image above the survey title)	Add new logo: Choose File No file chosen (Images wider than 600 pixels will be downsized to fit page.) If using a logo, hide survey title on survey page?			

• When you are finished editing your survey settings, scroll down to the bottom of the page and select "Save Changes." Your instrument is now enabled as a survey.

Send confirmation email (optional)? (Email the respondent when they complete the survey)	No V
	Save Changes
	Cancel